CE 0459

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Set for infusion of ophthalmic gas

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ARCEOLE • INDICATIONS : Endotamponade of the retina.

DESCRIPTION:

ARCEOLE is composed of a sterile unidose cylinder of 30ml non sterile gas for ophthalmic surgery and a sterile plastic graduated 50ml syringe mounted with a sterilizing 0.22µm filter and a device to connect it onto the cylinder. A 30G needle is supplied for the intraocular injection. A patient wrist band allows the identification of intraocular gas wearer. A set of 5 patient labels is provided to allow traceability.

The different accessories provided with ARCEOLE have been validated for the herunder detailed use, it is not advisable to substitute them by others.

COMPOSITION :

ARCEOLE is a sterile unidose cylinder of one of the following gases : SF₆ (Sulfur hexafluoride), C₂F₆ (Hexafluorohexane), C₃F₈ (Octafluoropropane).

- Each sterile unidose cylinder contains a pure non sterile gas in a sterile packaging.
- Nb : The gases included in ARCEOLE are non toxic, inert, non flammable, colorless and odorless compounds. Once implanted, they are not metabolized and gradually eliminated through bloodstream and the respiratory tract.
- Expansion and reabsorption rates based on bibliography data :

Gas	Expansion rate	Time allowed for expansion in days	Average tamponade duration in days	Non expansive concentration in %
SF ₆	1.9-2	2	10-15	20
C ₂ F ₆	3.3	3	30-35	16-17
C ₃ F ₈	4	3	55-65	12

PRECAUTIONS AND WARNINGS:

- Although non toxic those gases are competing with oxygen when inhaled. Avoid intentional inhalation.
- **It is imperative to stop the anaesthesia with nitrogen protoxide inhalation at least 15 mn before starting the gas injection into the patient eye.**
- The gases for vitreo retinal surgery are exclusively intended to be used by ophthalmologists trained to this technique.
- ARCEOLE is a medical device designed to inject an ophthalmic gas mixed with air, In case of injection of pure gas under the physician responsibility, The bibliography⁽¹⁾ recommends not to exceed the following quantity:

SF ₆	C ₂ F ₆	C ₃ F ₈
1.5 ml	1 ml	0.8 ml

- Avoid injection of pure gas in case of narrow irido-corneal angle, a pre-operative gonioscopy is recommended

⁽¹⁾JACOBS PM, TWOMEY JM, LEAVER PK. Behaviour of intraocular gases. Eye 1988, 2, 660-3

- After implantation, a treatment against ocular hypertension should be established and IOP daily monitored.
- Before and after injection of gas, monitoring of the retina central artery is required
- Perform regular follow up of glaucoma patients or patients with IOP, the elderly patients and those with diabetes.
- As a rule, it is necessary to check the ability of the patient to hold his head in the required position during the treatment.

COMPLICATIONS AND SIDE EFFECTS:

The complications described in the literature are those linked to vitreo-retinal surgery and particularly:

- sub-retinal passage, haemorrhage, haematoma, choroid tear, sub-conjunctival passage.
- Crystalline lens or retina touch with the needle during the injection.
- Risk of separation or capture of the gas bubble between the crystalline lens and the pars plana.
- Potential side effects are: hypertony, retinal tear, cataract, inflammation, endophtalmie, keratopathy and must be considered by the surgeon with regards to the expected benefits.

A PATIENT IMPLANTED WITH INTRAOCULAR GAS MUST ABSOLUTELY AVOID PRESSURE VARIATION (elevation, travel by plane, diving with or without isobar caisson, hyperbar caisson treatment)
A PATIENT IMPLANTED WITH INTRAOCULAR GAS MUST ABSOLUTELY AVOID NITROGEN PROTOXIDE ANAESTHESIA

In order to remind those side effects and warnings the patient is encouraged to wear the wrist band included in the package.

INSTRUCTIONS FOR USE:

WARNING the gas is not sterile and must be sterilized by passing through the sterilizing filter supplied with the product. Preparation must be performed extemporaneously.

- Position the syringe equipped with the filter and connector onto the protective seal of the valve and push vigorously the tip of the connector through the seal inside the valve.
- The gas is released inside the syringe and pushes back the plunger, it is sterilized by passing through the 0.22µm filter.
- Disconnect the syringe-filter part from the connector holding the cylinder and swinging the syringe with a right-left movement. (The filter and connector must stay connected to the syringe). Do not pull to disconnect.
- Adjust the plunger to the wished volume of pure gas, using the syringe graduation.

The table below provides help to evaluate the gas concentration:

Gas	Non expansive concentration in %	Volume of pure gas to fill a 50 ml syringe	Volume of air
SF ₆	20	10 ml	40 ml
C ₂ F ₆	16	8 ml	42 ml
C ₃ F ₈	12	6 ml	44 ml

- Take air in the syringe to obtain the desired proportion mix. The air is sterilized by passing through the syringe filter,
- Remove the filter from the syringe,
- Immediately mount a sterile injection canula with its protective cap, the mix is ready for intraocular injection.

STORAGE CONDITIONS:

Avoid high temperatures above 50°C/112°F

Store in dry, aerated area.

Do not use if the sterile pouch has been compromised.

Do not use after expiration date.

Single use

Do not re-sterilize

STERILIZATION:

NOTE: The gas contained in ARCEOLE is not sterile, only the cylinder and accessories are sterilized with ethylene oxide.

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